



Advancing Research with Combined Synchrotron Techniques

Lu Ma

QAS Lead Scientist

NSLS II, Brookhaven National Laboratory

03/12/2025



@BrookhavenLab

Outline

Introduction

- Overview of synchrotron X-ray characterization techniques
- Importance of combined techniques


Combining Techniques

- XAS + XRD: Comprehensive structural analysis
- XAS + DRIFTS: Catalyst insights


Pseudo grazing incidence measurement

- Total electron yield for surface analysis


Introduction




Multifunctional Additives for Next Generation Battery



X-ray Spectroscopy



X-ray Diffraction




Structure

Operando
Multi-modal
Synchrotron

Chemistry

Morphology

X-ray Microscopy




- Comprehensive Analysis
- Real-Time Monitoring
- Maximize data yield from valuable beamtime
- Reveal unexpected properties or phenomena
- Bridge gaps between research fields.

Introduction

Na-ion battery

Was that a slow reaction? Why did that happened?




During Charge:

Na⁺ extraction, oxidation at cathode

During Discharge:

Na⁺ insertion, reduction at cathode



Lack of Cross-Validation when using a single characterization technique!

Introduction




Coin cell in real life




Coin cell at the beamline

Modified with a hole in the middle to allow the X-ray penetration



76 pm

102 pm




J. Phys. Chem. Lett. 2015, 6, 11, 2081–2085

Combined XAS & XRD


X-ray diffraction (long-range order):

- Phase identification of crystalline materials
- Crystal structure
- Particle size, strain, and other microstructural properties





X-ray absorption spectroscopy (XAS) (short-range order):

- Electronic structure and oxidation state of specific elements
- Local coordination geometry and bond lengths





Combined XAS & XRD

NSLS-II
7-BM





MAX-IV Balder



J. Appl. Cryst. (2014). 47, 449–457.
J. Phys. Chem. C 2017, 121, 18202–18213.


Chem. Rec. 2019, 19, 1444–1456.
ACS Sustainable Chem. Eng. 2017, 5, 3631–3636.

Combined XAS & XRD









Combined XAS & XRD

- Case study 1a: Battery C/5 rate
XAS: Mn, Co, Ni K-edge
XRD: $\lambda = 0.9547 \text{ \AA}$




David Bock, BNL



ACS Appl. Mater. Interfaces 2021, 13, 43, 50920–50935.




Combined XAS & XRD

- Case study 1b: Fast- charging Battery



Tianyi Li, ANL


Lithium intensities after 900 fast charging cycles



Combined XAS & XRD

- Case study 1b: Fast-charging Battery
- Real time monitoring, 64X speed



6C rate, 2 spots
XAS: Ni K-edge
XRD: $\lambda = 0.8856 \text{ \AA}$



Unpublished

Combined XAS & XRD

- Case study 1b: Fast- charging Battery




Unpublished


Combined XAS & XRD

- Case study 1b: Fast- charging Battery

Bad spot of 1 layer
Presence of dead lithium






Good spot of 1 layer
No dead lithium





Combined XAS & XRD

- Case study 1b: Fast- charging Battery



Combined XAS & XRD

- Case study 2: Linkam stage heating





At 330°C (after temperature correction),
A layer structure begun to form, and Ni was reduced.


Under review

Combined XAS & XRD


- User's hydrothermal setup



Xin Zhang, PNNL



Combined XAS & DRIFTS



XAS provides information on the oxidation state, electronic structure, and local coordination environment around specific elements.


DRIFTS is sensitive to molecular vibrations and offers insight into surface species, functional groups, and adsorbed molecules on catalyst surfaces.

Combined XAS/DRIFTS: simultaneously monitor:



- Changes in the metal center (via XAS)
- Organic species or reactants/products (via DRIFTS) during a reaction.

Combined XAS & DRIFTS



ESRF ID24




APS 9-BM



QAS hutch C endstation -- DRIFTS



- Thermo-Nicolet iS-50 IR spectrometer and Harrick cell
- Plexiglass box to avoid influence from the air for DRIFTS measurement
- Transmission XAS only, fixed sample depth of 3mm
- Suitable for high-energy measurements but not for diluted samples or low-energy measurements below 12 keV.




Combined XAS & DRIFTS

- Science commissioning result

~5wt.% Ru

Lin Yuan,
Stanford University




Combined XAS & DRIFTS

- Science commissioning result

0.7wt%Fe samples


Reduced the depth to 1mm with a graphite insert

Dominik Wierzbicki,
TES, NSLS-II



Combined XAS & DRIFTS


- Fluorescence capability under development



Front: X-ray in





Back: Fluorescence




Other combined techniques with XAS


XRF mapping



Low  High
S concentration




XANES tomography



Pseudo grazing incidence -TEY



Surface sensitive measurement with TEY

In collaboration with TES
LBS Dr. Yonghua Du




With incident angle of 1 deg

- Surface nearest Cu-Cu1 bond distance increase
- Coordination number of the first two Cu-Cu bonds decrease





Pseudo grazing incidence -TEY

penetration depth for Cu ($\rho=8.96$) @ 9000.0eV





penetration depth for Cu ($\rho=8.96$) @ 9000.0eV



Surface Science 324 (1995) L371-L377

Pseudo grazing incidence -TEY

How water quickly effected the Cu surface and make the cutting much easier?



Take-home message

- **Combined synchrotron techniques in a single experiment at one beamline** provide complementary insights, offering a more comprehensive understanding of materials and their behaviors.
- However, integrating multiple techniques often requires **compromises**, as sample conditions need to accommodate both methods.
- **Careful experimental design** is crucial to ensure that conditions are feasible for both techniques.
- Despite the challenges, this approach significantly enhances **real-time, multi-dimensional characterization**.



Thank you!